



Model 9500/9550 Respiration/Heart Rate Monitor Technical Manual



**REQUITRON
MEDICAL INC**

MINNEAPOLIS, MINNESOTA 55447 U.S.A.

©Copyright 1991, Aequitron Medical,[®] Inc., 14800 28th Avenue North, Minneapolis, Minnesota 55447 U.S.A. All rights reserved. No portion of this manual may be copied, reproduced, or stored in any form without the express written permission of Aequitron Medical, Inc.

Aequitron Medical is a registered trademark of Aequitron Medical, Inc.

For more information: Contact your Aequitron representative for information on our full line of medical equipment and related services. Or, contact Aequitron directly. Phone: (800) 497-4979; (612) 557-9200. Fax: (612) 557-8200.

Contents

Chapter	Title	Page
1	Introduction	1-1
2	Description	2-1
	Front Panel	2-2
	Back Panel	2-3
	Accessories	2-4
3	Theory of Operation	3-1
	Power System	3-2
	Isolated Section	3-3
	Analog System	3-4
	Digital System	3-6
	9500/9550 System Block Diagrams	3-9
	9500/9550 Power Distribution Block Diagram	3-10
	ECG Channel Block Diagram	3-11
	Respiratory Channel/Lead Alarm Block Diagram	3-12
	Oximeter Block Diagram	3-13
	Microprocessor Hardware Block Diagram	3-14
	LEDs and Alarm Block Diagram	3-15
4	Inspection	4-1
	Equipment	4-2
	Inspection Procedure	4-3
	Infant Alarm Settings	4-6
	Adult Alarm Settings	4-7
	ECG Response	4-9
	Inverted Polarity	4-10
	Heart Rate = Resp.	4-11
	Lead Output	4-11
	Respiration Response	4-12
	Oxygen Desat	4-13
	Interface Test	4-13
	Generating a 9500/9550 Report Printout	4-14
	Summary Report Samples	4-15
	Charging Light	4-19
	Ground Wire Resistance	4-19
	Checksheet	
5	Disassembly	5-1

6	Schematics	6-1
	Final Assembly	6-2
	Digital Board	6-5
	Analog Board	6-9
	Analog Board (50Hz)	6-12
	Front Panel	6-15
	Front Panel Board	6-16
	Back Panel	6-18
	Back Panel Board	6-20
	Phone Jack Board	6-23
	Key Switch Board	6-25

Introduction

This manual presents technical and service information about Aequitron Medical Models 9500, 9550, and 9550/50 Respiration/Heart Rate monitors. This manual is intended for use by competent biomedical technicians who have been trained on this product by Aequitron Medical,[®] Inc. For information on the operation and application of the monitors, see the Dealer's, User's, and Reports Manuals.

Aequitron Medical believes the information herein is complete and accurate, but accepts no liability for errors, omissions, or misrepresentations.

AEQUITRON MEDICAL FURTHER DECLINES ANY WARRANTIES, EXPRESSED OR IMPLIED, FOR THE REPAIRED PRODUCT, INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

It is the user's responsibility to ensure that the product is properly maintained and that it is in safe and proper operating condition before it is put into use.

Refer any adjustments or procedures that exceed the scope of this manual to an Aequitron Technical Service Representative by calling:

(800) 497-3787

Chapters

Following is a list of this manual's chapters, with a summary of each. Use this list to guide you to the area of interest.

Theory of Operation A detailed review of the operation of the monitor's circuitry. This chapter describes the monitor down to the functional-block level.

Inspection A procedure to test and verify the operation of all major functions of the monitor. Use this procedure before and after any repairs to a monitor, and as part of periodic inspections of the monitors.

Disassembly A pictorial representation, with accompanying text, of the disassembly of a monitor. Disassembly is shown to the replaceable module level supported by Aequitron Medical.

Schematics Electrical schematics, component locators, and mechanical drawings for the major sub-assemblies of the monitors.

Conventions

Notes, Cautions, and Warnings mean the following throughout this technical manual:

Note Directions that make it easier to use or service the monitor.

Caution Directions that prevent damaging the equipment.

Warning Directions that warn of hazards to the patient, to a caregiver, or to service personnel.

Description

The monitor sounds a warning alarm and turns on the appropriate light on the front panel to indicate which alarm occurred. More than one alarm can occur at the same time. In this case, a light turns on for each alarm.

The monitor sounds alarms when it detects:

- A long pause in the patient's breathing effort
- A heart rate that is too fast or too slow
- An equipment problem

Your monitor is a warning system. It does not prevent problems with breathing or heart rate. It sounds an alarm and indicates where the problem is, but YOU must take the necessary action.


The front panel has:

- Lights for human and equipment alarms
- RESET button
- POWER switch
- Patient cable jack
- Green lights for breathing effort and heart beat
- Alarm speaker

The back panel has:

- Battery charger jack
- Green charger light
- Other jacks and controls for your doctor or homecare dealer

Front Panel



a Human Alarm Lights

Human alarms sound when the monitor detects a pause in the patient's breathing effort or a heart rate that is too fast or too slow.

During a human alarm the alarm beeps and a red light turns on to show you which alarm occurred.

b Equipment Alarm Lights

Equipment alarms sound when a problem occurs with the battery, with a connection to the patient, or with the monitor itself. During most equipment alarms the alarm sounds constantly. A yellow light turns on to show you which alarm occurred.

c Reset

Use this button to turn off the Human and LEAD alarm lights after an alarm condition is over. Also, hold this button in while moving the POWER switch to turn the monitor off.

d Power On/Off

To turn the monitor on, slide this switch to On. To turn the monitor off, hold in the RESET button and move the POWER switch to OFF. If you move the POWER switch to OFF without holding in the RESET button, an alarm will sound.

e Patient Cable Jack

Plug the patient cable into this jack.

f Respiration Light

This green light blinks once each time the monitor detects breathing effort.

g Heart Beat Light


This green light blinks once each time the monitor detects a heart beat.

h Alarm Speaker

The alarm sounds come from here.

WARNING Do not cover the alarm speaker (h), or you may not hear the alarm.

Back Panel



i Serial Number

j Battery Charger Jack

This is where you plug in the battery charger.

k Battery Charger Light

This green light turns on to indicate that the battery charger is properly connected.

l Interface

This output is used for recording respiration and heart rate. Other accessories may also be connected to this output.

m O₂ Sat/Rate (9550 series only)

On the Model 9550, this is the input for a pulse oximeter that your doctor may have you use with the monitor. Instructions for use are included with this device.

n ECG and RESP Outputs

These are outputs for accessories.

o Remote Alarm Jack

This jack is for an optional alarm.

p Controls



These controls are used with the display to set the alarm limits for breathing and heart rate. Never change these settings unless the doctor or homecare dealer tells you to.

q Function Display

When you first turn on your monitor, this display will show the limits for alarms and memory status. If the display shows the message, "Memory Full," follow your homecare dealer's instructions.

Accessories

These accessories pick up the patient's breathing and heart beat signals and carry them to your monitor. Use only Aequitron approved accessories.



Electrode Pads

The black electrode pads pick up the breathing and heart beat signals from the patient.


Rinse the flat sides of the pads with water each day. Do not use soap or detergent on the pads because both interfere with breath and heart beat signals.

Rinse all soap off the patient after a bath. Do not use oil, powder, or lotions on the patient where the electrodes will be placed. All these can interfere with picking up breath and heart beat signals.

Belt

The soft foam belt holds the electrode pads in place. The belt goes around the chest and is fastened with a Velcro strip.


Remove the electrode pads before washing the belt. You can wash belts by hand or in the washing machine using the delicate cycle. Use lukewarm water, a mild detergent, and no bleach. Rinse well to remove all detergent because leftover detergent can cause a skin rash. Air dry the belts. Never dry them in the automatic dryer because heat damages the inside of the belt.



Lead Wires

The lead wires pick up the breath and heart beat signals from the electrodes and carry them to the patient cable. Always hold the lead wires by the plastic ends when you are connecting or disconnecting them; otherwise, you can damage the wires inside the lead wires.

New lead wires come joined together. When you connect the patient to the monitor, pull gently to separate the lead wires as much as you need to.



Patient Cable

The patient cable picks up the breath and heart beat signals from the lead wires and carries the signals to the monitor. The cable has a locking connector where it attaches to the monitor. Hold the patient cable at the end when you are connecting or disconnecting it. Never pull on the cable's cord because you can damage the wires inside.

Theory of Operation

The monitors are intended for use with infants through adults. It detects and alarms for central apnea, tachycardia and bradycardia. Human events, equipment alarms and monitor on/off times are logged. Trend and waveform data for human events are also recorded. The 9550 series monitors also detect, record and validate oxygen desaturations when used in conjunction with an external pulse oximeter.

The monitor has an audible alarm (minimum of 79 dB(A) at 2 feet). It also has circuitry which detects the loss of any or all of the power rails. A watchdog circuit ensures that the microprocessor is running and is executing the software program correctly.

The monitor consists of five separate assemblies: the front panel, the analog board, the digital board, the back panel, and the battery. The phone jack board and key switch board are separate boards that plug into the back panel board.

For this theory, we will divide the monitor into four main parts: the power system, the isolated section, the analog system, and the digital system. The power system provides the power rails for the monitor. The isolated section provides patient isolation, supplies the drive current for respiration detection, and separates the ECG and respiratory signals. The analog section processes the ECG and respiratory signals, and digitizes them. The digital section performs calculations, detects events, and stores data.

Refer to the Block Diagram on page 3-9. On the block diagrams, the letters A, D, R, and F indicate the board where the particular circuit is located: A = analog board, D = digital board, R = rear (or back) panel, and F = front panel.

Power System

Refer to the Block Diagram on page 3-10.

The Model 9271 battery charger is a three-prong charger that provides 11.6 VDC at 180 mA. The charging (or operating) current is provided through diode CR1, which prevents a shorted input jack from damaging the device. The transistor, Z1, prevents voltage spikes (greater than ± 15 V) from damaging the unit or resetting the microprocessor.

The battery consists of eight nickel-cadmium (NiCad) cells connected in series. It is rated at 2000 mA·H (milli-amp hours). Fully charged battery voltage is between 11 and 12 VDC. The monitor will operate up to three days on a fully charged battery alone. The battery can be recharged in about 18 hours.

The battery voltage, V_B , passes through a parallel set of switches (Q11 and the front-panel POWER switch), after which it is known as V_A . V_A powers the alarm circuitry. Q11 will maintain power in the monitor (and set off the alarm) if the POWER switch is turned off without holding in the RESET switch. The fuse, F1, protects the rest of the circuitry from excessive current. The voltage on the protected side of the fuse is called V_P . V_M is derived from V_B and provides continuous power to the RAM chips and real time clock.

V_P provides the voltage for the various regulators in the monitor: +5.4 V for the Holter-type recorder, +5 V supplies, V_+ and V_- for the analog section. V_P is regulated down to V_S (+5 V) by VR2.

Note that circuit ground and chassis ground are tied together on the back panel board. They are also connected together (with back-to-back diodes) on the analog board.

Oscillator U4 is an astable multivibrator operating at approximately 30 kHz. The frequency is set by R34 and C19. The output at U4-10 has a 50% duty cycle. This signal drives Q1 which drives the primary of T1. T1 provides the negative voltage V_- , filtered by C30 and rectified by CR12. T1 also couples the 30 kHz through to the isolated section where it provides the power and carrier frequency for the patient signal.

Isolated Section

Power for the isolated section is derived from the 30 kHz provided by T1. The diodes, CR16-CR19, form a full-wave rectifier, which develops the positive and negative power for the isolated section. The isolated positive power rail is filtered by R71 and C41. The isolated negative power rail is filtered by R70 and C38.

Resistors R64 and R63 limit the 30 kHz current for the zener diodes (CR14 and CR15). These diodes set the drive voltage to approximately ± 10 VDC. C39 and C40 block any low frequency modulation of the patient drive. R65 and R66 limit the patient current to approximately 50 μ A. Since R65 and R66 are significantly larger than the patient impedance (typically 500 to 1000 ohms), the patient drive is essentially a constant-current source. Therefore, changes in the patient impedance amplitude-modulate the 30 kHz voltage across RA and LA. The respiration-modulated 30 kHz signal and the ECG appear across R67. The patient ECG signal can be thought of as a voltage source within the patient generating signals between 0.1 and 2 mV in amplitude. The respiration signal can be thought of as a small change in impedance (0.2 to 3 ohms) out of a total body impedance of 500 to 1000 ohms. The change in impedance amplitude-modulates the 30 kHz carrier signal.

R78, R79 and diodes CR23-CR30 provide input protection to the instrumentation amp. Surge arrestors E1, E2, and E3 provide protection to the circuitry if an electrostatic discharge is delivered to the patient or to the patient cable input.

The instrumentation amp U8 (pins 1,2,3, 5,6,7, 8,9,10) provides an extremely high input impedance, with a differential gain of 5 for both the ECG and respiratory signal.

C33 and R56 form a high pass filter. This couples the respiration-modulated 30 kHz signal to the primary of T2 and blocks the ECG signal. T2 couples the respiration signal to the non-isolated section.

C36 blocks DC while the network of L1, R74 and C37 provides a low-pass filter. This filters out the 30 kHz carrier leaving the ECG signal. U8-12,13,14 provides a gain of about nine, while C35 provides low-pass filtering with an f_c of 100 Hz.

Q2 and Q3 provide a full-wave chopper (at 30 kHz) that acts as the carrier for the ECG signal. T3 couples it to the non-isolated section of the analog board. R57 and R39 form a high impedance connection between isolated and non-isolated grounds to limit offset voltage build-up.

Analog System

The analog system consists of the ECG and respiration channels, and, in the 9550 series, the oximeter channel. The analog system processes the ECG and respiratory signals. It supplies signals for the monitor's outputs and the digital system's inputs.

ECG Channel

Refer to the Block Diagram on page 3-11.

R35 and C21 provide a 90° phase shift in the 30 kHz signal so it can be used to demodulate the ECG signal. The demodulator signal controls analog switch U5-3,4. The ECG signal is demodulated and filtered (R40 and C22), and DC blocked (C23).

U6-1,2,3 provides variable, non-inverting gain from 13 to 31, adjustable via R93. C42 rolls the signal off at about 50 Hz. R95 sets the DC offset for the rest of the ECG channel.

U6-12,13,14 and U6-5,6,7 form a line-frequency notch filter. U6-12,13,14 is a bandpass filter whose passband can be tuned to 50 or 60 Hz with R94 acting as the tuning element. The output of U6-12,13,14 is 180° out-of-phase with the 50 or 60 Hz noise signal at U6-1. U6-5,6,7 sums the signals from U6-1 and TP9 (the inverted signal at 50 or 60 Hz) to provide a signal at TP10 that has limited AC line-frequency noise. R92 is used to adjust the summing of the two signals.

U6-8,9,10 R56, C44, R55, C28, R52 and C43 form a third-order Butterworth low pass filter with unity gain. The corner frequency is about 40 Hz. The output of this section is sent to the back panel Interface connector (J5) through R48 (ECG High). The signal is also divided by R47 and R46 and again sent to the Interface connector J5 (ECG Low) for use by Holter-type recorders.

U7-1,2,3 acts as a bandpass filter (center frequency of about 40 Hz) with gain of about 24. This filter is lightly damped (i.e., the output will ring) and will detect heartbeats whether RA/LA leads are correct or reversed.

U7-8,9,10 is a peak detector. Positive voltages (heartbeats) at U7-8 will pass through CR22 and be held on the capacitor C48. The voltages held on C48 decay off through R96 if there is no positive voltage at U7-10.

U7-12,13,14 and V- provide the ECG signal with a DC offset. R23 and C5 filter the ECG signal at U3-2 such that it can be compared with the unfiltered version at U3-3. The output, U3-1, changes from low to high when a heartbeat is detected.

The ECG signal (high or low) arrives on the digital board at connector J12-5 and controls the gate of Q1, which is fed into Port A. The microprocessor receives the signal via its data bus.

LEAD Alarm Refer to the Block Diagram on page 3-12.

U1-5,6,7 acts as a peak detector providing a voltage at TP2 that corresponds to the impedance across the patient leads. This voltage is divided by R2 and R104, then sent to input 3 of the A/D via J12-6. It's then transmitted over the data bus to the microprocessor. Nominally, patient impedances greater than 1800 ohms will cause a lead alarm.

Respiration Channel Refer to the Block Diagram on page 3-12. R35 and C21 provide a 90° phase shift in the 30 kHz signal. This is used to demodulate the respiratory signal. The demodulator signal controls the analog switch U5-1,2. R28 and C14 filter the respiratory signal. U2-1,2,3 provides non-inverting gain of gain of about 18.

R11, C13, R15, C16, R29, C15 and U2-5,6,7 form a third-order Butterworth low-pass filter with unity gain. The corner frequency is about 5 Hz.

U2-12,13,14 provides an inverting gain of between 300 and 600, depending on the adjustment of R20.

CR1 and CR2 are 5.1 V zener diodes. The current flowing through these two diodes and R7 keeps the voltage output of U2-14 from going into saturation, hence, allowing it to recover quickly from large swings in its output. If the output of U2-7 (between C17 and R19) makes large swings, the diodes CR3 and CR4, the zeners CR1 and CR2, and the signal inversion in U2-12,13,14 combine to keep the output (U2-14) from saturating.

U2-8,9,10 functions as a bandpass filter and as a buffer for the fast recovery amp (U2-12,13,14). C12 and R9 give a lower corner frequency of 0.07 Hz and R14 and C7 give an upper corner frequency of about 9 Hz. The output of U2-8 goes out to the back panel Interface connector, J5, via R50 as the RESP High signal. The output also is filtered by the C24, C25, R30, R32 network and again goes out to the Interface connector, J5 (Resp Filtered). In addition, the output

is divided by R49 and R51, DC blocked by C29 and delivered to the back panel as the RESP Low signal for use by Holter-type recorders.

R1 and R8 provide a DC offset at U1-14 and at U1-8. C8 and C9 roll off the respiratory signal at 10 Hz. The signal is filtered by R27 and C11, and buffered by U1-1,2,3. From here the signal goes to input 0 of the A/D (U10) via P12-4. It is converted and transferred into the microprocessor via the data bus.

Oximeter Channel Refer to the Block Diagram on page 3-13.

The oximeter sends out two signals. One is a voltage proportional to the percentage of oxygen saturation (O₂ SAT). The second is a voltage proportional to the heart rate (Pulse Rate). These two signals are fed directly back out on the Interface connector (J5).

U2-1,2,3 and U2-12,13,14 act as buffer op amps with nominal gains of 4 to prepare the Pulse Rate and O₂ SAT signal for the A/D. O₂ SAT is on input 4 of the A/D and Pulse Rate goes in on input 5. The signals are converted and transferred to the microprocessor via the data bus.

Digital System

Refer to the Block Diagram on page 3-14.

The digital system consists of a microprocessor with its support chips, the system memory, a watchdog circuit, a real-time clock, a serial port, and input/output hardware. The monitor provides battery back-up for portions of the digital system. The digital system performs calculations, detects events, and stores data in the monitor's memory.

Microprocessor The 65C02 microprocessor is a low-power device using an 8 bit data bus and 16 bit address bus. It runs at a 1 MHz clock rate provided by crystal Y1 and the U4 oscillator circuit.

U14 provides time and date information to the microprocessor via the data bus.

U10 is an eight input, eight bit A/D converter. A0, A1, A2 are decoded to select one of the eight inputs. It uses the microprocessor 1 MHz clock as its timing reference. It uses a 5 V reference which means that eight bits can define 256 different voltage levels leading to each distinct voltage level being about 19 mV.

U2 is a Versatile Interface Adaptor (VIA). This chip provides two bidirectional, eight bit ports for the microprocessor. The VIA also provides the timing for the interrupt cycle.

Two ports, U7 and U8, are output ports used for controlling the LEDs and the LCD via the data bus.

U9 provides the interface to communicate with a personal computer via a serial communications port in a standard RS-232C format. In general, the six signals sent and received by this chip allow communications between the microprocessor and the computer.

Memory The two RAM chips, U11 and U12, provide storage space for events that are recorded by the monitor. About 44 kB of these chips' 64 kB of available storage space is used.

The software for the system is stored in the socketed EPROM chip, U13. The serial number for the unit is also stored in the EPROM so that the checksum for each EPROM is unique.

Watchdog The watchdog circuit causes an alarm if the microprocessor fails to operate properly. The watchdog circuit consists of two retriggerable one-shots (U17), and a flip-flop (U16-4,5,6,8,9,10). On start-up, C32 keeps U16-5 low and U16-4 high for about four seconds. The high at U16-4 turns on Q12. Q12 turns on Q11 and ensures that power will be supplied to the circuitry (V_B to V_A) even if the on/off switch is turned off. R36 and C31 provide a delay at start-up to ensure that the MONITOR LED is turned off initially. Every 16 mSec (nominally) the microprocessor toggles the watchdog circuit input through Port B (PB-7). This signal triggers the first one-shot output (U17-7) to go low for 12 mSec. The output of the first one-shot (U17-7) triggers the second one-shot output (U17-10) to go high for 50 mSec. The continuous high output of the second one-shot (U17-10) is sent to the active low SET line (U16-8) of the flip-flop. If the microprocessor toggles the watchdog input faster than 12 mSec or slower than 50 mSec, then the one-shots time out, U16-8 goes low, U16-10 goes high (which ultimately turns on the MONITOR LED on the front panel) and the audible alarm is turned on via U16-11,12,13 and Q10.

LED and Alarm Drivers

Refer to the Block Diagram on page 3-15. The **RESPIRATION LED** blinks each time the monitor detects breathing effort. The microprocessor generates a signal on the data bus. It is sent out through Port B (PB-3). A high turns on Q6. Q6 drives the **RESPIRATION LED** and the Resp sense line of the Interface connector (J5).

The **HEARTBEAT LED** blinks each time the monitor detects a heart beat. The microprocessor generates a signal on the data bus. It is sent through Port B (PB-4). A high turns on Q5. Q5 drives the **HEART BEAT LED** and the ECG sense line of the Interface connector (J5).

Human alarm conditions (apnea, heart slow, or heart fast) cause a loud alarm at one second intervals. Equipment failure conditions (battery, lead or monitor) cause a loud, continuous alarm. The soft audio signal notifies the user of actions taken (for example, when a button is pressed), of a battery caution alarm, or of an invalid limits alarm.

The audio alarm can be turned on by the watchdog circuit. Normally, U16-4 is a high which means U16-11 is low (assume U16-12 is high). Q10 is off and no alarm sounds. If the watchdog times out, U16-4 goes low, U16-11 goes high, Q10 turns on and drives the alarm.

The loud alarm can be turned on by a low from the microprocessor, sent out through Port B (PB-5), putting a low at U16-12 and turning on Q10. Q10 drives the audio alarm.


The soft audio alarm is turned on by a low generated by the microprocessor, sent out through Port B (PB-6). This is inverted by U23-1,2, which turns on Q9 with a high at its gate. Q9 drives the audio alarm through R45, which limits the current flow through the alarm, hence, reducing the alarms sound output level.

A loud audio alarm will also be generated if Vp or Vs are pulled low. If this happens, U15-7 goes high, turning on Q8, which drives the audio alarm. This will occur if the internal fuse (F1) blows.


A fault sensed by the watchdog timer will turn on the **MONITOR LED** as described earlier. The final drive circuitry (U5-13,11 and U23-12,13) is also driven by a signal generated by the microprocessor and fed through Port B (PB-2). The **MONITOR LED** will flash for an invalid limits alarm.

Alarm conditions are detected through software in the monitor. The drive for the alarm LEDs is provided by the microprocessor data bus through Port C. A low at the output pin of Port C will turn on the front panel LED.


95xx System Block Diagram




95XX Power Distribution




ECG Channel




Respiratory Channel / Lead Alarm






Microprocessor Hardware



LEDs and Alarm



Inspection

The following inspection procedure tests all major functions of 9500 or 9550 series monitors. Perform this procedure before and after you repair any monitor. Standard operating procedure may also require that you perform this test at additional intervals to ensure continued, optimum operation of the monitor.

Throughout this procedure, enter the test results called for on the Inspection Data Sheet. You will find a sample of this sheet at the end of this chapter. You may copy this sample as needed for use during all monitor checkout procedures.

Perform the procedure before repairing a monitor, entering your results in the Initial Test column of the data sheet. This may help you isolate the defective assembly. After repair, repeat the procedure, entering your results in the Final Test column of the data sheet.

Headings in the procedure correspond to the tests listed on the data sheet. Perform all the tests, in the order listed. If a monitor fails any test, or any part of any test, that monitor fails the entire procedure.

Warning

Do not use any monitor that does not pass all parts of the inspection procedure. For assistance, contact Aequitron at (800) 497-3787.

If the MONITOR alarm LED comes on at any time, except during the power-up LED test, reject the unit.

Equipment

To complete this procedure, you will need the following equipment:

- Aequitron Model 8310 or Model 9310 Simulator
- Aequitron Model 8237 Patient Cable
- Aequitron Model 8245 Lead Wires
- Aequitron Model 9271 Battery Charger
- Aequitron Model 8200-14 70 dB Cap
- Final Test Junction Box¹
- Co-axial cable
- Digital multimeter (DMM)
- Oscilloscope (scope) with digital storage—
The oscilloscope must have a time-base range of at least 50 msec. per division to 0.5 sec. per division, and a vertical input range of at least 5 mV per division to 2 V per division.

Caution

An Aequitron Model 8310/50 or Model 9310/50 simulator must be used to test a 9550/50 monitor. A 50 Hz power supply should be used. An Aequitron Model 9271/50 battery charger is to be used with a 9550/50 monitor.

Note The monitors are designed to meet AAMI "Safe Current Limits for Electromedical Apparatus," July 9, 1985.

¹ The Final Test Junction Box is available only to persons who complete the Model 9500/9550 repair training.

Inspection Procedure

Serial # Read the serial number on the back of the monitor and enter it on the Inspection Data Sheet. Write your name and today's date in the appropriate blanks on the data sheet.


Case Assembly Visually inspect the monitor for proper case assembly. All parts of the case must fit securely, with no gaps between components. Wires must not protrude between case components.

Mechanical Verify the operation of all switches and the ball. Inspect all hardware, including jacks. Switches must operate smoothly throughout their ranges. The ball must secure in an upright position. All hardware must be tight and secure. All jacks must accept and hold their respective plugs without undue resistance.

Alarm With no connections to the monitor, and with the face of the audible alarm free of any obstructions, turn the monitor on. The alarm must sound with a clear, audible tone. Turn the monitor off.

Figure 4.1

Test equipment connections
 Beginning with the Power up test, connect the equipment as shown. Note that the patient cable is not yet connected to the monitor.



Power Up

Connect the equipment as shown in Figure 4.1. Note that the patient cable is not yet connected to the monitor. You may place a Model 8200-14 70 dB cap over the alarm face to muffle the audible alarm. **You must remove the cap before you return the monitor to service.** Set the equipment controls as listed in Table 4.1.

Turn on the monitor. All alarm LEDs (front panel) must light and turn off sequentially. There will be an approximately two second delay before the audible lead alarm sounds. Hold down the Remote Alarm button on the Junction Box. The Junction Box alarm must sound while you press the button. Release the Audible Alarm button. The monitor's audible alarm and **LEAD** LED remain on.

Connect the patient cable to the monitor. The audible alarm turns off. Then the monitor alarms for apnea and heart slow within 12 seconds. The audible alarm beeps once each second. The **HEART BEAT** and **RESPIRATION** LEDs are off; the **APNEA**, **HEART SLOW**, and **LEAD** LEDs are on.

On the simulator, press Resp. Rate to end the apnea, and set ECG Amplitude to 1 mV. The monitor's alarm must be silent. The **HEART BEAT** and **RESPIRATION** LEDs blink. The **APNEA**, **HEART SLOW**, and **LEAD** LEDs remain on. On the monitor, press **RESET**. The **APNEA**, **HEART SLOW**, and **LEAD** LEDs turn off. Turn off the monitor.

Table 4.1 Initial test equipment settings
Set the equipment controls as listed before the power up test.

Equipment	Settings
Simulator	Power: On; Respiration Rate: Apnea (Inf.); Respiration Impedance: 1 Ω ; Electrode Impedance: 500 Ω ; ECG Rate: 165 BPM; ECG Ampl: 0 mV
Monitor	Power: Off; Apnea: 10 sec; Heart Fast: Off; Heart Slow: 40 BPM;
Junction Box	Selector: Recorder Power

Recorder Power Set the junction box selector to Recorder. With the scope/DMM, measure the recorder power output. The voltage level must be 5.4 ± 0.5 VDC.

Note: If monitor is being checked for annual recertification, the battery should be replaced.

2 YRS

Battery Output Change the selector on the Junction Box to Battery. With the scope, measure the battery output. The battery voltage must be greater than 10.0 VDC.

Resp. Output Set the Junction Box selector to Resp. With the scope, verify the output is 2.0 ± 0.2 V_{p-p}, with a maximum offset of 1 VDC.

ECG Output Set the Junction Box selector to ECG. With the scope, verify the output is 3.0 ± 0.3 V_{pk}, from the ECG baseline to the top of the R-wave.

Infant Alarm Settings Set the monitor as follows in Table 4.2:

Table 4.2 Infant Alarm Settings
Set the monitor and simulator.

Equipment	Settings		
Monitor	New Patient		
	Infant		
	Monitor Alarm Limits	Apnea	10 sec
		Heart Fast	300 BPM
Heart Slow		130 BPM	
Simulator	ECG rate		165 BPM
	ECG amplitude		1 mV
	Resp rate		30 BPM
	Resp impedance		1.0 Ω
	Electrode impedance		500 Ω

Table 4.3 Verify Infant Alarms
Verify that all alarms are working correctly at each setting.

Set:	Setting	Verify Alarm:	Action
ECG rate	315 BPM	Heart Fast alarm activates	none
	290 BPM	Audible alarm is off and HEART FAST LED is on.	Press RESET to turn LED off.
	115 BPM	Heart Slow alarm activates.	none
	140 BPM	Audible alarm is off and HEART SLOW LED is on.	Press RESET to turn LED off.
Simulator	Inf Apnea, and press Enter.	Apnea alarm after 8.5 to 11.5 seconds.	none
	Press Resp Rate on simulator.	Audible alarm is off and APNEA LED is on.	Press RESET to turn LED off.

Adult Alarm Settings

Table 4.4 Adult Alarm Settings
Set the monitor and simulator.

Equipment	Settings		
Monitor	New Patient		
	Adult		
	Patient Name	Enter your initials	
	Patient Number	1234	
	Monitor Alarm Limits	Apnea	60 sec
		Heart Fast	100 BPM
		Heart Slow	30 BPM
	Event Log Limits	Apnea	45 sec
		Heart Fast	100 BPM
		Heart Slow	30 BPM
		O2 Sat (9550 series only)	80%
	Displayed Events (9550 series only)	Priority 1	Apnea
		Priority 2	Desat
		Priority 3	HR Slow
Verify that "SAT %" = $83 \pm 2\%$ and "Pulse" = 165 ± 3 BPM. (9550 series only) Pulse can be adjusted with Pulse knob on the junction box.			
Time/Date	current time and date		
Simulator	ECG rate	95 BPM	
	ECG amplitude	1 mV	
	Resp rate	10 BPM	
	Resp impedance	1.0 Ω	
	Electrode impedance	500 Ω	

Adult Alarm Settings, continued

Table 4.5 **Verify Adult Alarms**
Verify that all alarms are working correctly at each setting.

Set:	Setting	Verify Alarm:	Action
ECG rate	105 BPM	Heart Fast alarm activates	none
	95 BPM	Audible alarm is off and HEART FAST LED is on.	Press RESET to turn LED off.
	25 BPM	Heart Slow alarm activates.	none
	35 BPM	Audible alarm is off and HEART SLOW LED is on.	Press RESET to turn LED off.
Simulator	Inf Apnea, and press Enter.	Apnea alarm after 58.5 to 61.5 seconds. (Start this timing after the last Resp pulse.)	none
	Press Resp Rate on simulator.	Audible alarm is off and APNEA LED is on.	Press RESET to turn LED off.

ECG Response

Table 4.6

ECG Response Settings
 Make the following changes to the monitor settings:

Equipment	Settings		
Monitor	Monitor Alarm Limits	Apnea	10 sec
		Heart Fast	off
		Heart Slow	30 BPM
	Event Log Limits	Apnea	10 sec
		Heart Fast	off
		Heart Slow	30 BPM
		O ₂ Sat (9550 series only)	80%
Junction Box Switch	ECG		

Table 4.7

ECG Response Verification
 Set the simulator for the following ECG rates and amplitudes.

Equipment	ECG Rate	ECG Amplitude	Verify:
Simulator	35 BPM	0.2 mV	Verify 10 beats at each setting. Each setting should have no missing beats or double beats. Verify that no alarms are activated.
	315 BPM	0.2 mV	
	315 BPM	5.0 mV	
	35 BPM	5.0 mV	

Inverted Polarity

Table 4.8 **Inverted Polarity Connections**
Make the following changes to the simulator lead connections:

Equipment	Leads	Verify
Simulator	Disconnect the RA and LA leads from the simulator.	Lead alarm acitvates
	Connect the RA lead to the LA simulator socket. Connect the LA lead to the RA simulator socket.	Audible alarm is off. Press RESET to turn LEAD LED off.

Table 4.9 **Inverted Polarity Verification**
Set the simulator for the following ECG rates and amplitudes.

Equipment	ECG Rate	ECG Amplitude	Verify:
Simulator	35 BPM	5 mV	Verify 10 beats at each setting. Each setting should have no missing beats or double beats. Verify that no alarms are activated.
	315 BPM	5 mV	
	315 BPM	0.5 mV	
	35 BPM	0.5 mV	
	35 BPM	0.2 mV	
	215 BPM	0.2 mV	
Lead Connections	Return the RA lead to the RA simulator socket. Return the LA lead to the LA simulator socket.		

Heart rate = Resp.

Table 4.10 Heart rate = Resp.
Check the monitor and simulator by following the steps in this table.

Equipment	Settings:	Verify
1. Simulator	ECG Rate = 95 BPM	
	Resp Rate = 75 BPM	
	Resp Impedance = 0.5 Ω	
2. Monitor	Press RESET .	
3. Simulator	ECG Rate = 75 BPM	The monitor's APNEA LED and audible alarm must turn on.
4. Simulator	Resp Rate = 50 BPM	The monitor's audible alarm must be silent, The APNEA LED must remain on.
5. Monitor	Press RESET .	The APNEA LED must turn off.

Lead Output

Table 4.11 Lead Output
Set the simulator for the following Electrode Impedance settings and verify the alarms.

Electrode Impedance Setting:	Verify:	Action
3000 Ω	Lead Alarm activates	
1500 Ω	Audible alarm turns off.	Press RESET to turn off the LEAD LED
1000 Ω	Normal operation of the monitor (no alarms).	
500 Ω	Normal operation of the monitor (no alarms).	

Respiration Response

Table 4.14 **Respiration Response**
Set the junction box and simulator for the following initial settings.

Equipment	Settings:
Junction Box switch	Resp
Simulator	ECG rate = 165 BPM
	ECG Amplitude = 1.0 mV
	Electrode Impedance = 500 Ω

Table 4.15 **Respiration Response**
Set the simulator for the following Respiration rates and Respiration Impedance settings.

Equipment	Resp Rate	Resp Impedance	Verify
Simulator	50 BPM	0.2 Ω	Verify 6 breaths at each setting. Each setting should have no missing breaths or double breaths. Verify that no alarms are activated.
	20 BPM	0.2 Ω	
	15 BPM	0.4 Ω	
	10 BPM	0.5 Ω	
	75 BPM	0.5 Ω	
	75 BPM	5.0 Ω	
	10 BPM	5.0 Ω	
Simulator	Resp Rate = 30 BPM	Resp Impedance = 1 Ω	

Oxygen DeSat

Table 4.15 Oxygen DeSat

Equipment	Setting:
Junction Box switch (9550 series only)	Press the O ₂ Sat switch for approximately 10 seconds, then release.

Interface Test

Transfer the data from the monitor to a personal computer that has the Report Generator software.

Page 4-14 outlines the basics of transferring data from the monitor to a computer. Refer to the Report Generator Manual for further details.

Verify that events simulated during testing were accurately transferred by reviewing the printed Summary Report and Event Log.

Examples of a 9550 Summary Report and Event Log are on page 4-15 and 4-16.

Examples of a 9500 Summary Report and Event Log are on page 4-17 and 4-18.

Generating a 9500/9550 Report/Printout

Computer	9500/9550 Monitor
1. 25 pin D connector cable between 9500/9550 and computer	
2. Turn on computer, and type CD\RG . Press Enter .	
3. Turn on printer.	
4. Insert formatted patient diskette.	
5. Type MAIN .	
	6. Turn on while holding "Function" and "Up" buttons.
7. Type R for "Receive Monitor Data."	
8. Type D for "Direct Connect."	
	9. Press "Down" button.
10. Press the Space Bar.	
	11. Three beeps means it's o.k.
12. Type E for "Examine Patient Data."	
13. Press Enter .	
14. Type S to view a summary report.	
15. Type P to print summary report.	
16. Type E to view an event log.	
17. Type P to print an event log.	
18. Type L to leave the program.	

Summary Report 9550 Sample

```

PRELIMINARY REPORT      ** MODEL 9550 **      SUMMARY REPORT

MONITOR INFORMATION      YOUR INITIALS      Adult
                          1234

** EVENT RECORDING TIME **      Time      Date

Event Recording Start      16:14:44      06/30/93
Event Recording Stop      16:28:10      06/30/93
Total Recording Time      14 Mins

Total Monitor On Time      12 Mins
% of Total Recording Time  85 %

** APNEAS **

Monitor Alarm Apneas >= 10 seconds      2
Event Logging Apneas >= 10 seconds      2

Longest Apnea lasted 69 secs at 16:17:03 on 06/30
Min Heart Rate (3 beat Avg) : 34 BPM      Min O2 Sat : Invalid

** BRADYCARDIAS **

Monitor Alarm Bradycardias (<= 30 BPM      1
Event Logging Bradycardias (<= 30 BPM      1

Lowest Heart Rate lasted 14 secs at 16:16:40 on 06/30
Min Heart Rate (3 beat Avg) : 23 BPM      Min O2 Sat : Invalid

** TACHYCARDIAS **

Monitor Alarm Tachycardias      OFF
Event Logging Tachycardias      OFF

Highest Heart Rate lasted 7 secs at 16:16:24 on 06/30
Max Heart Rate (3 beat Avg) : 105 BPM      Min O2 Sat : Invalid

** DESATURATIONS **

Event Logging Desaturations      1

Lowest O2 Saturation lasted 13 secs at 16:25:50 on 06/30
Min Heart Rate (3 beat Avg) : 165 BPM      Min O2 Sat : 77 %

** EQUIPMENT ALARMS **

Loose Lead Alarms      3 Alarms
Low Battery Alarms      -

** DISPLAYED EVENTS **      )) Are Not Reviewed ((

Priority #1 (20) Apnea      2 Found
Priority #2 (10) Desat      1 Found
Priority #3 (10) Brady      1 Found

Copyright (c) 1989,92 Aquitron Medical Inc. Rel 040-029-00433075
    
```

Event Log 9550 Sample

```

AQUITRON MEDICAL INC.      ** MODEL 9500 **      EVENTS LOG
Patient Name: YOUR INITIALS      Number: 1234

Event Type      Time      Date      Duration      HR (BPM)
                (3 beat Avg)
1. Switch Change Monitor Alarm      Event Logger
   Apnea Limit      60 secs      45 secs
   Tachy Limit      100 BPM      100 BPM
   Brady Limit      30 BPM      30 BPM

2. Monitor On      16:17:59      07/31      2 Mins
3. Tachycardia      16:18:36      07/31      8 Secs      105      *
4. Bradycardia      16:18:51 ++ 07/31      14 Secs      23      *
5. Apnea            16:19:24 ++ 07/31      1 Min 12 Secs 35      *
6. Monitor Off      16:20:32      07/31      1 Min
7. Switch Change Monitor Alarm      Event Logger
   Apnea Limit      10 secs      10 secs
   Tachy Limit      OFF          OFF
   Brady Limit      30 BPM      30 BPM

8. Monitor On      16:21:46      07/31      6 Mins
9. Loose Lead      16:23:13      07/31      4 Secs
10. Loose Lead      16:25:08      07/31      4 Secs
11. Apnea           16:25:35 ++ 07/31      18 Secs      75      *
12. Loose Lead      16:26:20      07/31      4 Secs
13. Monitor Off      16:27:59      07/31      11 Mins

End of Logged Events      Copyright (c) 1989,92 Aquitron Medical Inc.
    
```

Summary Report 9500 Sample

PRELIMINARY REPORT ** MODEL 9500 ** SUMMARY REPORT

MONITOR INFORMATION YOUR INITIALS Adult
1234

** EVENT RECORDING TIME ** Time Date
Event Recording Start 16:17:59 07/31/93
Event Recording Stop 16:27:59 07/31/93
Total Recording Time 10 Mins

Total Monitor On Time 9 Mins
% of Total Recording Time 90 %

**** APNEAS ****

Monitor Alarm Apneas >= 10 seconds 2
Event Logging Apneas >= 10 seconds 2

Longest Apnea lasted 72 secs at 16:19:24 on 07/31
Min Heart Rate (3 beat Avg) : 35 BPM

**** BRADYCARDIAS ****

Monitor Alarm Bradycardias <= 30 BPM 1
Event Logging Bradycardias <= 30 BPM 1

Lowest Heart Rate lasted 14 secs at 16:18:51 on 07/31
Min Heart Rate (3 beat Avg) : 23 BPM

**** TACHYCARDIAS ****

Monitor Alarm Tachycardias OFF
Event Logging Tachycardias OFF

Highest Heart Rate lasted 8 secs at 16:18:36 on 07/31
Max Heart Rate (3 beat Avg) : 105 BPM

**** EQUIPMENT ALARMS ****

Loose Lead Alarms 3 Alarms
Low Battery Alarms -

**** DISPLAYED EVENTS ****)) Are Not Reviewed ((

Priority #1 (20) Apnea 2 Found
Priority #2 (10) Brady 1 Found
Priority #3 (10) Brady -

Copyright (c) 1989,92 Requitron Medical Inc. Rel 040-029-00462322

*REV SOFTWARE COMPUTER
REV II MONITOR
MONITOR SOURCE*

Event Log 9500 Sample

AEQUITRON MEDICAL INC. ** MODEL 9550 ** EVENTS LOG

Patient Name: YOUR INITIALS Number: 1234

Event Type	Time	Date	Duration	HR (BPM) (3 beat Avg)	O2 (%)	
1. Switch Change	Monitor Alarm		Event Logger			
Apnea Limit	60 secs		45 secs			
Tachy Limit	100 BPM		100 BPM			
Brady Limit	30 BPM		30 BPM		O2 Sat Limit 80%	
2. Monitor On	16:14:44	06/30	3 Mins			
3. Tachycardia	16:16:24	06/30	7 Secs	105	INV	*
4. Bradycardia	16:16:40 ++	06/30	14 Secs	23	INV	*
5. Apnea	16:17:02 ++	06/30	1 Min 9 Secs	34	INV	*
6. Monitor Off	16:18:22	06/30	1 Min			
7. Switch Change	Monitor Alarm		Event Logger			
Apnea Limit	10 secs		10 secs			
Tachy Limit	OFF		OFF			
Brady Limit	30 BPM		30 BPM		O2 Sat Limit 80%	
8. Monitor On	16:19:14	06/30	8 Mins			
9. Loose Lead	16:20:58	06/30	4 Secs			
10. Loose Lead	16:22:18	06/30	3 Secs			
11. Apnea	16:22:49 ++	06/30	15 Secs	75	INV	*
12. Loose Lead	16:23:28	06/30	3 Secs			
13. Desaturation	16:25:50 ++	06/30	13 Secs	165	77	
14. Monitor Off	16:28:10	06/30	6 Mins			

End of Logged Events Copyright (c) 1989,92 Aequitron Medical Inc.

Charging Light

Turn off the monitor and disconnect all equipment. Plug battery charger into monitor (back panel) and verify that the green charge LED lights. Disconnect the battery charger.

Ground Wire Resistance

Using a digital multimeter, measure the resistance between the shell of the Interface "D" connector on the back panel and the exposed metal of either the ECG or Resp output on the back panel. Resistance should be one ohm or less.

Warning

Do not use the exposed metal of the REMOTE ALARM connector on the back panel for the groundwire resistance check.

Sub-D Conn. = tighten
if fails

Model _____

Model 9500/9550 Inspection Data Sheet

Serial No. _____

Initial Inspection Pass Fail

Final Inspection Pass Fail

Inspected By: _____

Inspected By: _____

Inspection Date: _____

Inspection Date: _____

Test	Initial Test		Result	Final Test		Result	Spec.
	Pass/Fail			Pass/Fail			
Serial #	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Case assembly	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Mechanical	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Alarm	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Power up	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Recorder power	<input type="checkbox"/>	<input type="checkbox"/>	VDC	<input type="checkbox"/>	<input type="checkbox"/>	VDC	5.4 ± 0.5 VDC
Battery output	<input type="checkbox"/>	<input type="checkbox"/>	VDC	<input type="checkbox"/>	<input type="checkbox"/>	VDC	>10.0 VDC
Resp. output	<input type="checkbox"/>	<input type="checkbox"/>	Vp-P	<input type="checkbox"/>	<input type="checkbox"/>	Vp-P	2.0 ± 0.2 Vp-P
ECG output	<input type="checkbox"/>	<input type="checkbox"/>	Vpk	<input type="checkbox"/>	<input type="checkbox"/>	Vpk	3.0 ± 0.3 Vpk
Infant Alarm Settings	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Adult Alarm Settings	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
ECG Response	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Inverted Polarity	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Heart Rate = Resp.	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Lead Output	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Respiration Response	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Oxygen Desat	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		9550 Only
Interface Test	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Charging light	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		
Ground Wire Resistance	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		≤ 1 ohm

Notes:

Disassembly

The following is a description of the disassembly process for a Model 9500 or Model 9550 monitor. It breaks the monitor into sub-assemblies that you can replace with Aequitron-supplied replacement parts. Disassembly beyond the level shown here will void all warranties and support contracts for that monitor.

Caution

You must always disassemble monitors only at a properly grounded workstation. Always wear a ground strap when disassembling monitors.

You will need the following tools to disassemble the monitor:

- A large flat-blade screwdriver
- A small knife
- One #2 cross-blade screwdriver
- One #0 cross-blade screwdriver
- One 1/4" nut driver
- One 3/16" nut driver

Figure 5.1 Protective Panel
 With a flat-blade screwdriver, remove the nylon screw from the back of the monitor top. Swing down the top of the protective panel. Pull the protective panel off of the monitor's back panel.




Figure 5.2 Feet and Ball
 Turn the monitor upside-down. With a small knife, remove the silicon plugs from the feet. Use a #2 cross-blade screwdriver to remove the screw from each foot. Remove the feet and the ball from the bottom of the monitor.





Figure 5.3

Remove the bottom and set it beside the unit. Disconnect the battery cable from the rear panel. Disconnect the four ground wires from the two studs on the battery cover. Lift off the two side panels.

**Figure 5.4**

Remove seven screws from the analog board. Disconnect J14, J15, and J19 from the analog board. Lift out the analog board.




Figure 5.5

Remove seven studs and lift out the shield. Remove the next seven studs and lift out the digital board with the front and rear panels.





Figure 5.6

Disconnect J13 from the digital board and remove the front panel. Lift the rear panel to disconnect J8 from the digital board.



Schematics

Following are the mechanical drawings, component locators, and electrical schematics for the final assembly of the monitors, and for their major subassemblies. For information on the disassembly of the monitors, see the Disassembly Chapter of this manual. For information on replacement parts, contact Aequitron's Technical Service Department by calling:


(800) 497-3787

The following service kits are used to replace defective assemblies in the monitors:

Monitor	Service Kit	Board
9500	1061	Digital Board
	1057	Analog Board
	1059	Front Panel
9550	1061	Digital Board
	1057	Analog Board
	1060	Front Panel
9550/50	1061	Digital Board
	1058	Analog Board
	1060	Front Panel

Note Back panels are not field-serviceable. Units with defective back panels must be returned to Aequitron Medical for service.






Final Assembly




INTERCONNECT DIAGRAM

- ⚠️ INSERT PROGRAMMED MODULE ON DIGITAL BOARD AT LOCATION UTIL.
 - ⚠️ PLACE INSULATOR ITEM 38 BETWEEN BATTERY PACK, ITEM 21, AND CASE BOTTOM, ITEM 2.
 - ⚠️ CENTER LABEL, ITEM 16, BETWEEN MOUNTING FEET AS SHOWN.
 - ⚠️ EARTH GROUND CONNECTED THROUGH STANDOFF.
 - ⚠️ GROUND STRAPS.
 - 3. RESISTANCE BETWEEN BATTERY CASE AND SHIELD (ITEM 32) MUST BE LESS THAN 1 OHM.
 - ⚠️ COMETER CONNECTOR TO APPEAR ON 02-003210-001 ONLY, REF.
1. NEWS AND SECTIONS SIMPLIFIED FOR CLARITY ON THIS DRAWING.
- NOTES:


← protective panel
p/w








Digital Board

PART USED	REVISION
CS1	CR9, CR10, CR11
CR17	A1
CR18	CR12
CR19	CR13
CR20	CR14
CR21	CR15
CR22	CR16
CR23	CR17
CR24	CR18
CR25	CR19
CR26	CR20
CR27	CR21
CR28	CR22
CR29	CR23
CR30	CR24
CR31	CR25
CR32	CR26
CR33	CR27
CR34	CR28
CR35	CR29
CR36	CR30
CR37	CR31
CR38	CR32
CR39	CR33
CR40	CR34
CR41	CR35
CR42	CR36
CR43	CR37
CR44	CR38
CR45	CR39
CR46	CR40
CR47	CR41
CR48	CR42
CR49	CR43
CR50	CR44
CR51	CR45
CR52	CR46
CR53	CR47
CR54	CR48
CR55	CR49
CR56	CR50
CR57	CR51
CR58	CR52
CR59	CR53
CR60	CR54
CR61	CR55
CR62	CR56
CR63	CR57
CR64	CR58
CR65	CR59
CR66	CR60
CR67	CR61
CR68	CR62
CR69	CR63
CR70	CR64
CR71	CR65
CR72	CR66
CR73	CR67
CR74	CR68
CR75	CR69
CR76	CR70
CR77	CR71
CR78	CR72
CR79	CR73
CR80	CR74
CR81	CR75
CR82	CR76
CR83	CR77
CR84	CR78
CR85	CR79
CR86	CR80
CR87	CR81
CR88	CR82
CR89	CR83
CR90	CR84
CR91	CR85
CR92	CR86
CR93	CR87
CR94	CR88
CR95	CR89
CR96	CR90
CR97	CR91
CR98	CR92
CR99	CR93
CR100	CR94
CR101	CR95
CR102	CR96
CR103	CR97
CR104	CR98
CR105	CR99
CR106	CR100
CR107	CR101
CR108	CR102
CR109	CR103
CR110	CR104
CR111	CR105
CR112	CR106
CR113	CR107
CR114	CR108
CR115	CR109
CR116	CR110
CR117	CR111
CR118	CR112
CR119	CR113
CR120	CR114
CR121	CR115
CR122	CR116
CR123	CR117
CR124	CR118
CR125	CR119
CR126	CR120
CR127	CR121
CR128	CR122
CR129	CR123
CR130	CR124
CR131	CR125
CR132	CR126
CR133	CR127
CR134	CR128
CR135	CR129
CR136	CR130
CR137	CR131
CR138	CR132
CR139	CR133
CR140	CR134
CR141	CR135
CR142	CR136
CR143	CR137
CR144	CR138
CR145	CR139
CR146	CR140
CR147	CR141
CR148	CR142
CR149	CR143
CR150	CR144
CR151	CR145
CR152	CR146
CR153	CR147
CR154	CR148
CR155	CR149
CR156	CR150
CR157	CR151
CR158	CR152
CR159	CR153
CR160	CR154
CR161	CR155
CR162	CR156
CR163	CR157
CR164	CR158
CR165	CR159
CR166	CR160
CR167	CR161
CR168	CR162
CR169	CR163
CR170	CR164
CR171	CR165
CR172	CR166
CR173	CR167
CR174	CR168
CR175	CR169
CR176	CR170
CR177	CR171
CR178	CR172
CR179	CR173
CR180	CR174
CR181	CR175
CR182	CR176
CR183	CR177
CR184	CR178
CR185	CR179
CR186	CR180
CR187	CR181
CR188	CR182
CR189	CR183
CR190	CR184
CR191	CR185
CR192	CR186
CR193	CR187
CR194	CR188
CR195	CR189
CR196	CR190
CR197	CR191
CR198	CR192
CR199	CR193
CR200	CR194
CR201	CR195
CR202	CR196
CR203	CR197
CR204	CR198
CR205	CR199
CR206	CR200
CR207	CR201
CR208	CR202
CR209	CR203
CR210	CR204
CR211	CR205
CR212	CR206
CR213	CR207
CR214	CR208
CR215	CR209
CR216	CR210
CR217	CR211
CR218	CR212
CR219	CR213
CR220	CR214
CR221	CR215
CR222	CR216
CR223	CR217
CR224	CR218
CR225	CR219
CR226	CR220
CR227	CR221
CR228	CR222
CR229	CR223
CR230	CR224
CR231	CR225
CR232	CR226
CR233	CR227
CR234	CR228
CR235	CR229
CR236	CR230
CR237	CR231
CR238	CR232
CR239	CR233
CR240	CR234
CR241	CR235
CR242	CR236
CR243	CR237
CR244	CR238
CR245	CR239
CR246	CR240
CR247	CR241
CR248	CR242
CR249	CR243
CR250	CR244
CR251	CR245
CR252	CR246
CR253	CR247
CR254	CR248
CR255	CR249
CR256	CR250
CR257	CR251
CR258	CR252
CR259	CR253
CR260	CR254
CR261	CR255
CR262	CR256
CR263	CR257
CR264	CR258
CR265	CR259
CR266	CR260
CR267	CR261
CR268	CR262
CR269	CR263
CR270	CR264
CR271	CR265
CR272	CR266
CR273	CR267
CR274	CR268
CR275	CR269
CR276	CR270
CR277	CR271
CR278	CR272
CR279	CR273
CR280	CR274
CR281	CR275
CR282	CR276
CR283	CR277
CR284	CR278
CR285	CR279
CR286	CR280
CR287	CR281
CR288	CR282
CR289	CR283
CR290	CR284
CR291	CR285
CR292	CR286
CR293	CR287
CR294	CR288
CR295	CR289
CR296	CR290
CR297	CR291
CR298	CR292
CR299	CR293
CR300	CR294
CR301	CR295
CR302	CR296
CR303	CR297
CR304	CR298
CR305	CR299
CR306	CR300
CR307	CR301
CR308	CR302
CR309	CR303
CR310	CR304
CR311	CR305
CR312	CR306
CR313	CR307
CR314	CR308
CR315	CR309
CR316	CR310
CR317	CR311
CR318	CR312
CR319	CR313
CR320	CR314
CR321	CR315
CR322	CR316
CR323	CR317
CR324	CR318
CR325	CR319
CR326	CR320
CR327	CR321
CR328	CR322
CR329	CR323
CR330	CR324
CR331	CR325
CR332	CR326
CR333	CR327
CR334	CR328
CR335	CR329
CR336	CR330
CR337	CR331
CR338	CR332
CR339	CR333
CR340	CR334
CR341	CR335
CR342	CR336
CR343	CR337
CR344	CR338
CR345	CR339
CR346	CR340
CR347	CR341
CR348	CR342
CR349	CR343
CR350	CR344
CR351	CR345
CR352	CR346
CR353	CR347
CR354	CR348
CR355	CR349
CR356	CR350
CR357	CR351
CR358	CR352
CR359	CR353
CR360	CR354
CR361	CR355
CR362	CR356
CR363	CR357
CR364	CR358
CR365	CR359
CR366	CR360
CR367	CR361
CR368	CR362
CR369	CR363
CR370	CR364
CR371	CR365
CR372	CR366
CR373	CR367
CR374	CR368
CR375	CR369
CR376	CR370
CR377	CR371
CR378	CR372
CR379	CR373
CR380	CR374
CR381	CR375
CR382	CR376
CR383	CR377
CR384	CR378
CR385	CR379
CR386	CR380
CR387	CR381
CR388	CR382
CR389	CR383
CR390	CR384
CR391	CR385
CR392	CR386
CR393	CR387
CR394	CR388
CR395	CR389
CR396	CR390
CR397	CR391
CR398	CR392
CR399	CR393
CR400	CR394
CR401	CR395
CR402	CR396
CR403	CR397
CR404	CR398
CR405	CR399
CR406	CR400
CR407	CR401
CR408	CR402
CR409	CR403
CR410	CR404
CR411	CR405
CR412	CR406
CR413	CR407
CR414	CR408
CR415	CR409
CR416	CR410
CR417	CR411
CR418	CR412
CR419	CR413
CR420	CR414
CR421	CR415
CR422	CR416
CR423	CR417
CR424	CR418
CR425	CR419
CR426	CR420
CR427	CR421
CR428	CR422
CR429	CR423
CR430	CR424
CR431	CR425
CR432	CR426
CR433	CR427
CR434	CR428
CR435	CR429
CR436	CR430
CR437	CR431
CR438	CR432
CR439	CR433
CR440	CR434
CR441	CR435
CR442	CR436
CR443	CR437
CR444	CR438
CR445	CR439
CR446	CR440
CR447	CR441
CR448	CR442
CR449	CR443
CR450	CR444
CR451	CR445
CR452	CR446
CR453	CR447
CR454	CR448
CR455	CR449
CR456	CR450
CR457	CR451
CR458	CR452
CR459	CR453
CR460	CR454
CR461	CR455
CR462	CR456
CR463	CR457
CR464	CR458
CR465	CR459
CR466	CR460
CR467	CR461
CR468	CR462
CR469	CR463
CR470	CR464
CR471	CR465
CR472	CR466
CR473	




SPARE GATES






Analog Board





△ TABS ON T1 ARE OPTIONAL.
 1. UNLESS OTHERWISE SPECIFIED:
 ALL DIODES ARE 1N4148
 ALL RESISTOR VALUES ARE IN OHMS,
 ALL 5% RESISTORS ARE 1/4W
 ALL CAPACITOR VALUES ARE IN
 MICROFARADS.


NOTES:

PARTS USED (OTHER ITEMS)	
Q1	1N4148
Q2	1N4148
Q3	1N4148
Q4	1N4148
Q5	1N4148
Q6	1N4148
Q7	1N4148
Q8	1N4148
Q9	1N4148
Q10	1N4148
Q11	1N4148
Q12	1N4148
Q13	1N4148
Q14	1N4148
Q15	1N4148
Q16	1N4148
Q17	1N4148
Q18	1N4148
Q19	1N4148
Q20	1N4148
Q21	1N4148
Q22	1N4148
Q23	1N4148
Q24	1N4148
Q25	1N4148
Q26	1N4148
Q27	1N4148
Q28	1N4148
Q29	1N4148
Q30	1N4148
Q31	1N4148
Q32	1N4148
Q33	1N4148
Q34	1N4148
Q35	1N4148
Q36	1N4148
Q37	1N4148
Q38	1N4148
Q39	1N4148
Q40	1N4148
Q41	1N4148
Q42	1N4148
Q43	1N4148
Q44	1N4148
Q45	1N4148
Q46	1N4148
Q47	1N4148
Q48	1N4148
Q49	1N4148
Q50	1N4148
Q51	1N4148
Q52	1N4148
Q53	1N4148
Q54	1N4148
Q55	1N4148
Q56	1N4148
Q57	1N4148
Q58	1N4148
Q59	1N4148
Q60	1N4148
Q61	1N4148
Q62	1N4148
Q63	1N4148
Q64	1N4148
Q65	1N4148
Q66	1N4148
Q67	1N4148
Q68	1N4148
Q69	1N4148
Q70	1N4148
Q71	1N4148
Q72	1N4148
Q73	1N4148
Q74	1N4148
Q75	1N4148
Q76	1N4148
Q77	1N4148
Q78	1N4148
Q79	1N4148
Q80	1N4148
Q81	1N4148
Q82	1N4148
Q83	1N4148
Q84	1N4148
Q85	1N4148
Q86	1N4148
Q87	1N4148
Q88	1N4148
Q89	1N4148
Q90	1N4148
Q91	1N4148
Q92	1N4148
Q93	1N4148
Q94	1N4148
Q95	1N4148
Q96	1N4148
Q97	1N4148
Q98	1N4148
Q99	1N4148
Q100	1N4148





Analog Board (50Hz)



NOTES:


△ TABS ON T1 ARE OPTIONAL.

1. UNLESS OTHERWISE SPECIFIED:
 ALL RESISTOR VALUES ARE IN OHMS,
 1/4W, 1%
 ALL CAPACITORS ARE 1/4W
 ALL CAPACITOR VALUES ARE IN
 MICROFARADS.

LAST USED DELETED ITEMS	
R104	
C104	
C105	
C106	
C107	
C108	
C109	
C110	
C111	
C112	
C113	
C114	
C115	
C116	
C117	
C118	
C119	
C120	
C121	
C122	
C123	
C124	
C125	
C126	
C127	
C128	
C129	
C130	
C131	
C132	
C133	
C134	
C135	
C136	
C137	
C138	
C139	
C140	
C141	
C142	
C143	
C144	
C145	
C146	
C147	
C148	
C149	
C150	
C151	
C152	
C153	
C154	
C155	
C156	
C157	
C158	
C159	
C160	
C161	
C162	
C163	
C164	
C165	
C166	
C167	
C168	
C169	
C170	
C171	
C172	
C173	
C174	
C175	
C176	
C177	
C178	
C179	
C180	
C181	
C182	
C183	
C184	
C185	
C186	
C187	
C188	
C189	
C190	
C191	
C192	
C193	
C194	
C195	
C196	
C197	
C198	
C199	
C200	
C201	
C202	
C203	
C204	
C205	
C206	
C207	
C208	
C209	
C210	
C211	
C212	
C213	
C214	
C215	
C216	
C217	
C218	
C219	
C220	
C221	
C222	
C223	
C224	
C225	
C226	
C227	
C228	
C229	
C230	
C231	
C232	
C233	
C234	
C235	
C236	
C237	
C238	
C239	
C240	
C241	
C242	
C243	
C244	
C245	
C246	
C247	
C248	
C249	
C250	
C251	
C252	
C253	
C254	
C255	
C256	
C257	
C258	
C259	
C260	
C261	
C262	
C263	
C264	
C265	
C266	
C267	
C268	
C269	
C270	
C271	
C272	
C273	
C274	
C275	
C276	
C277	
C278	
C279	
C280	
C281	
C282	
C283	
C284	
C285	
C286	
C287	
C288	
C289	
C290	
C291	
C292	
C293	
C294	
C295	
C296	
C297	
C298	
C299	
C300	
C301	
C302	
C303	
C304	
C305	
C306	
C307	
C308	
C309	
C310	
C311	
C312	
C313	
C314	
C315	
C316	
C317	
C318	
C319	
C320	
C321	
C322	
C323	
C324	
C325	
C326	
C327	
C328	
C329	
C330	
C331	
C332	
C333	
C334	
C335	
C336	
C337	
C338	
C339	
C340	
C341	
C342	
C343	
C344	
C345	
C346	
C347	
C348	
C349	
C350	
C351	
C352	
C353	
C354	
C355	
C356	
C357	
C358	
C359	
C360	
C361	
C362	
C363	
C364	
C365	
C366	
C367	
C368	
C369	
C370	
C371	
C372	
C373	
C374	
C375	
C376	
C377	
C378	
C379	
C380	
C381	
C382	
C383	
C384	
C385	
C386	
C387	
C388	
C389	
C390	
C391	
C392	
C393	
C394	
C395	
C396	
C397	
C398	
C399	
C400	
C401	
C402	
C403	
C404	
C405	
C406	
C407	
C408	
C409	
C410	
C411	
C412	
C413	
C414	
C415	
C416	
C417	
C418	
C419	
C420	
C421	
C422	
C423	
C424	
C425	
C426	
C427	
C428	
C429	
C430	
C431	
C432	
C433	
C434	
C435	
C436	
C437	
C438	
C439	
C440	
C441	
C442	
C443	
C444	
C445	
C446	
C447	
C448	
C449	
C450	
C451	
C452	
C453	
C454	
C455	
C456	
C457	
C458	
C459	
C460	
C461	
C462	
C463	
C464	
C465	
C466	
C467	
C468	
C469	
C470	
C471	
C472	
C473	
C474	
C475	
C476	
C477	
C478	
C479	
C480	
C481	
C482	
C483	
C484	
C485	
C486	
C487	
C488	
C489	
C490	
C491	
C492	
C493	
C494	
C495	
C496	
C497	
C498	
C499	
C500	
C501	
C502	
C503	
C504	
C505	
C506	
C507	
C508	
C509	
C510	
C511	
C512	
C513	
C514	
C515	
C516	
C517	
C518	
C519	
C520	
C521	
C522	
C523	
C524	
C525	
C526	
C527	
C528	
C529	
C530	
C531	
C532	
C533	
C534	
C535	
C536	
C537	
C538	
C539	
C540	
C541	
C542	
C543	
C544	
C545	
C546	
C547	
C548	
C549	
C550	
C551	
C552	
C553	
C554	
C555	
C556	
C557	
C558	
C559	
C560	
C561	
C562	
C563	
C564	
C565	
C566	
C567	
C568	
C569	
C570	
C571	
C572	
C573	
C574	
C575	
C576	
C577	
C578	
C579	
C580	
C581	
C582	
C583	
C584	
C585	
C586	
C587	
C588	
C589	
C590	
C591	
C592	
C593	
C594	
C595	
C596	
C597	
C598	
C599	
C600	
C601	
C602	
C603	
C604	
C605	
C606	
C607	
C608	
C609	
C610	
C611	
C612	
C613	
C614	
C615	
C616	
C617	
C618	
C619	
C620	
C621	
C622	
C623	
C624	
C625	
C626	
C627	
C628	
C629	
C630	
C631	
C632	
C633	
C634	
C635	
C636	
C637	
C638	
C639	
C640	
C641	
C642	
C643	
C644	
C645	
C646	
C647	
C648	
C649	
C650	
C651	
C652	
C653	
C654	
C655	
C656	
C657	
C658	
C659	
C660	
C661	
C662	
C663	
C664	
C665	
C666	
C667	
C668	
C669	
C670	
C671	
C672	
C673	
C674	
C675	
C676	
C677	
C678	
C679	
C680	
C681	
C682	
C683	
C684	
C685	
C686	
C687	
C688	
C689	
C690	
C691	
C692	
C693	
C694	
C695	
C696	
C697	
C698	
C699	
C700	
C701	
C702	
C703	
C704	
C705	
C706	
C707	
C708	
C709	
C710	
C711	
C712	
C713	
C714	
C715	
C716	
C717	
C718	
C719	
C720	
C721	
C722	
C723	
C724	
C725	
C726	
C727	
C728	
C729	
C730	
C731	
C732	
C733	
C734	
C735	
C736	
C737	
C738	
C739	
C740	
C741	
C742	
C743	
C744	
C745	
C746	
C747	
C748	
C749	
C750	
C751	
C752	
C753	
C754	
C755	
C756	
C757	
C758	
C759	
C760	
C761	
C762	
C763	
C764	
C765	
C766	
C767	
C768	
C769	
C770	
C771	
C772	
C773	
C774	
C775	
C776	
C777	
C778	
C779	
C780	
C781	
C782	
C783	
C784	
C785	
C786	
C787	
C788	
C789	
C790	
C791	
C792	
C793	
C794	
C795	
C796	


Analog Board (50Hz)

Model 9500/9550: Technical Manual





0.55 MAX
COMPONENT HT.
(TRANSFORMER)
T.P. E KEAMPT

Front Panel



Front Panel Board



LAST USED: C5
CR10
J18
Q1
R11
SW2



DELETED: J1-15

3. ALL PARTS TO BE SEATED SQUARELY ON PCB. MAX HEIGHT FOR Q1 IS .375.


2. MASK FOR J17 ONLY.

1. UNLESS OTHERWISE SPECIFIED, ALL RESISTOR VALUES ARE IN OHMS, 1/8W, 1%; ALL CAPACITOR VALUES ARE IN MICROFARADS.


NOTES:




DIGITAL BOARD




Back Panel



- ▲ 1. INSERT KEYING PLUG, ITEM 17, INTO D-CONNECTOR, PINS 1 AND 24. BREAK OFF TAB.
 - ▲ 2. CONNECTOR, ITEM 3, TO BE SOLDERED AND DEFLUXED AFTER ITEMS 1 AND 2 ARE FASTENED TOGETHER (REF).
 - ▲ 3. POSITION #4 EYELET END OF GROUND STRAP, ITEM 13, PARALLEL TO EDGE OF PANEL AS SHOWN.
 - ▲ 4. THIS FEATURE TO APPEAR ON -001 ONLY.
- NOTES:



VIEW A-A



Back Panel Board


LAST USED:	DELETED:	REF.:
C20	CR10, 11	J3
CR12	P1-7	J5
R10	Q1	C19
Q1	R15, 16	C2
R17		
V1		
V2		
Z1		


- ▲ PLACE FERRITE SUBSTRATE OVER PINS OF CONNECTORS AFTER SOLDERING. TACK IN PLACE WITH RTV AFTER SOLDERING THE D-CONNECTOR.
- ▲ THESE PINS MUST BE SEALED SQUARELY WITH PINS IN A PERPENDICULAR ORIENTATION WITH RESPECT TO PCB (SEE FIG. 3), J3, J5-7, P8, J9, J10, C20.
- ▲ MAIN ON COMPONENT MOUNT FOR COMPONENTS SOLDERED TO ITEM 1 IS .500 INCHES. J3, CR2 AND J5 ARE SEPARATELY SPECIFIED ON SHEET 3. THIS BOARD IS NOT TO BE REDESIGNED WITHOUT INCLUDED WITH THIS HEAVY RESTRICTION.


UNLESS OTHERWISE SPECIFIED: ALL CAPACITORS ARE IN OHMS, ALL RESISTORS ARE IN OHMS, ALL CAPACITOR VALUES ARE IN MICROFARADS.

NOTES:


LCD DISPLAY








Phone Jack Board




2. ITEM 2, EARTH TERMINAL CONNECTION, TO BE ATTACHED TO REAR PANEL STUD.

1. THESE ITEMS MUST BE SEATED SQUARELY ON PCB (ITEM 1). SKEWING OR CANTING OF PARTS IS NOT ACCEPTABLE.

NOTES:



Key Switch Board




△ THESE ITEM MUST BE SEATED SQUARELY ON PCB (ITEM 1). P6 & P7 PINS SHOULD BE PERPENDICULAR TO PCB SURFACE. S2 - S5 CANNOT BE TILTED.

1. UNLESS OTHERWISE SPECIFIED, ALL RESISTOR VALUES ARE IN OHMS, 1/8W, 1%; ALL CAPACITOR VALUES ARE IN MICROFARADS.

NOTES:

Key Switch Board

Model 9500/9550: Technical Manual



Model _____

Serial No. _____

Model 9500/9550 Inspection Data Sheet


Initial Inspection Pass Fail Final Inspection Pass Fail

Inspected By: _____ Inspected By: _____

Inspection Date: _____ Inspection Date: _____

Test	Initial Test		Result	Final Test		Spec.
	Pass/Fail	Result		Pass/Fail	Result	
Serial #	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Case assembly	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Mechanical	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Alarm	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Power up	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Recorder power	<input type="checkbox"/>	<input type="checkbox"/>	VDC	<input type="checkbox"/>	<input type="checkbox"/>	VDC 5.4 ± 0.5 VDC
Battery output	<input type="checkbox"/>	<input type="checkbox"/>	VDC	<input type="checkbox"/>	<input type="checkbox"/>	VDC >10.0 VDC
Resp. output	<input type="checkbox"/>	<input type="checkbox"/>	V _{p-p}	<input type="checkbox"/>	<input type="checkbox"/>	V _{p-p} 2.0 ± 0.2 V _{p-p}
ECG output	<input type="checkbox"/>	<input type="checkbox"/>	V _{pk}	<input type="checkbox"/>	<input type="checkbox"/>	V _{pk} 3.0 ± 0.3 V _{pk}
Infant Alarm Settings	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Adult Alarm Settings	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
ECG Response	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Inverted Polarity	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Heart Rate = Resp.	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Lead Output	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Respiration Response	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Oxygen Desat	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	9550 Only
Interface Test	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Charging light	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	
Ground Wire Resistance	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	≤ 1 ohm

Notes:



MODEL 9500/9550 MONITOR ALARM QUESTION & ANSWER SUMMARY

Please refer to pages 11, 23 and 29 of the 9500 Dealers Manual and page 22 of the 9500 Users Manual for further discussion of these issues.

Q: What type of Monitor alarms are there?

A: 2 types. In the first, the monitor light (on the front panel) is on continuously and the audible alarm is sounding continuously. In the second, the monitor light is blinking once per second and the audible alarm is sounding once every 30 seconds.

Q: Under what conditions will the first condition occur (Monitor light on continuously, audible alarm sounding continuously)?

A: The microprocessor in the monitor continuously checks itself (via circuitry called a 'watchdog'). If the circuitry detects a potential error condition it will turn on the monitor light/alarm. Sometimes these conditions are transient and at others times the damage is permanent. The following are conditions under which this type of alarm might occur: electronics parts failure, battery discharged to extremely low levels, a large electrostatic discharge or a violent drop (of the monitor).

Q: What should I do when a continuous Monitor alarm occurs?

A: Turn the monitor off and then on again. If the problem was transient then the monitor will go through a normal start-up sequence and start to pick up heart and respiration signals. If the monitor goes into a second continuous Monitor alarm condition, we recommend that the monitor be returned for further evaluation. If you have any doubts at all about the integrity of the monitor please return it for evaluation.

Q: Under what circumstances will the second condition occur (the monitor light is blinking (once per second) and the audible alarm is sounding once every 30 seconds)?

A: The monitor stores its settings in several different locations in its internal memory. It regularly checks to see that these various locations agree. If the locations don't agree then it returns all the settings to their 'default' values (specifically, infant default settings) and it turns on the blinking Monitor light and beeps once every 30 seconds. The monitor continues to operate properly, however it may not be using the prescribed settings. A disruption to part (or all) of the internal memory could cause this to happen. The disruptions could be a result of any of the following: electronics parts failure (intermittent), battery discharged to extremely low levels, electrostatic discharge (of large magnitude) or a violent drop (of the monitor).

Q: What should I do if the monitor goes into this default condition (the monitor light is blinking once per second and the audible alarm is sounding once every 30 seconds)?

A: Under these conditions, the monitor is quite capable of continuing to safely monitor the patient, however, you must be aware that the alarm and logging limit settings may not be the prescribed ones. We recommend that you 'clear' the memory and re-enter the patient name/number and settings. Again, if you have any doubts as to whether or not the monitor is operating correctly, please call one of our tech service reps or return the monitor for further evaluation.

(1)

The first part of the document discusses the importance of maintaining accurate records of all transactions and activities. It emphasizes that this is essential for ensuring transparency and accountability in the organization's operations.

Furthermore, it highlights the need for regular audits and reviews to identify any discrepancies or areas for improvement. This process should be conducted in a systematic and thorough manner to ensure the highest level of accuracy.

In addition, the document stresses the importance of clear communication and collaboration between all departments. This will help to ensure that everyone is working towards the same goals and objectives, and that any issues are identified and resolved promptly.

Overall, the document provides a comprehensive overview of the organization's current state and outlines the key areas for focus and improvement. It is hoped that these recommendations will be implemented effectively to ensure the organization's long-term success and growth.

(2)

The second part of the document focuses on the financial aspects of the organization. It provides a detailed analysis of the current budget and identifies areas where costs can be reduced without compromising the quality of services provided.

It also discusses the importance of maintaining a strong financial position and the need for careful planning and forecasting. This will help to ensure that the organization has sufficient resources to meet its obligations and invest in future growth opportunities.

Furthermore, the document highlights the need for transparency in financial reporting and the importance of providing accurate and timely information to stakeholders. This will help to build trust and confidence in the organization's financial management.

In conclusion, the document provides a clear and concise summary of the organization's financial situation and offers practical recommendations for improvement. It is hoped that these measures will be implemented to ensure the organization's financial stability and long-term success.

(3)

The final part of the document discusses the human resources aspect of the organization. It emphasizes the importance of attracting and retaining top talent and provides recommendations for developing a strong and motivated workforce.

Overall, the document provides a comprehensive overview of the organization's current state and outlines the key areas for focus and improvement. It is hoped that these recommendations will be implemented effectively to ensure the organization's long-term success and growth.